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
Tare riferimento a:

Health Protection Agency publications:

"Laboratory reporting to the HPA. A guide for diagnostic laboratories".


"Hospital infection control : Guidance on the control of infection in hospitals".

Local guidelines including Infection Control Policy and Memorandum of Understanding.

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Bibliografia

1. European Parliament. UK Standards for Microbiology Investigations (SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes".
2. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices. 7-12-1998. p. 1-37.
3. Coia JE, Duckworth GJ, Edwards Di, Farrington M, Fry C, Humphreys H, et al. Guidelines for the control and prevention of methicillin-resistant *Staphylococcus aureus* (MRSA) in healthcare facilities. *J Hosp Infect* 2006;63 Suppl 1:S1-44.
4. Scottish Infections Standards and Strategies (SISS) Group. Guidance for the management of methicillin-resistant *Staphylococcus aureus*. The Royal College of Physicians of Edinburgh and the Royal College of Physicians and Surgeons of Glasgow. 2006.
5. Department of Health. Saving Lives: a delivery programme to reduce Healthcare Associated Infection including MRSA. Screening for Methicillin-resistant *Staphylococcus aureus* (MRSA) colonisation: A strategy for NHS trusts: a summary of best practice. Department of Health. 2006.
6. Jevons P. "Celbenin" - resistant Staphylococci. *British Medical Journal* 1961;1:124.
7. Schaeffler S, Jones D, Perry W, Ruvinskaya L, Baradet T, Mayr E, et al. Emergence of gentamicin- and methicillin-resistant *Staphylococcus aureus* strains in New York City hospitals. *J Clin Microbiol* 1981;13:754-9.
8. Pavillard R, Harvey K, Douglas D, Hewstone A, Andrew J, Collopy B, et al. Epidemic of hospital-acquired infection due to methicillin-resistant *Staphylococcus aureus* in major Victorian hospitals. *Med J Aust* 1982;1:451-4.
9. Oliveira DC, Tomasz A, de Lencastre H. The evolution of pandemic clones of methicillin-resistant *Staphylococcus aureus*: identification of two ancestral genetic backgrounds and the associated *mec* elements. *Microb Drug Resist* 2001;7:349-61.
10. Cookson BD. Nosocomial antimicrobial resistance surveillance. *J Hosp Infect* 1999;43 Suppl:S97-103.
11. Speller DC, Johnson AP, James D, Marples RR, Charlett A, George RC. Resistance to methicillin and other antibiotics in isolates of *Staphylococcus aureus* from blood and cerebrospinal fluid, England and Wales, 1989-95. *Lancet* 1997;350:323-5.
12. CDSC. *Staphylococcus aureus* bacteraemia: England and Wales, 2001. *CDR Weekly* 2002;12:1-17.
13. Johnson AP, Davies J, Guy R, Abernethy J, Sheridan E, Pearson A, et al. Mandatory surveillance of methicillin-resistant *Staphylococcus aureus* (MRSA) bacteraemia in England: the first 10 years. *J Antimicrob Chemother* 2012;67:802-9.
14. Enoch DA, Cargill JS, Sismey A, Karas JA. MRSA surveillance in a UK district hospital: measuring clinical isolates with MRSA is more useful than measuring MRSA bacteraemias. *J Hosp Infect* 2011;79:287-91.

15. Jernigan JA, Clemence MA, Stott GA, Titus MG, Alexander CH, Palumbo CM, et al. Control of methicillin-resistant *Staphylococcus aureus* at a university hospital: one decade later. *Infect Control Hosp Epidemiol* 1995;16:686-96.
16. Stamm AM, Long MN, Belcher B. Higher overall nosocomial infection rate because of increased attack rate of methicillin-resistant *Staphylococcus aureus*. *Am J Infect Control* 1993;21:70-4.
17. David MZ, Daum RS. Community-associated methicillin-resistant *Staphylococcus aureus*: epidemiology and clinical consequences of an emerging epidemic. *Clin Microbiol Rev* 2010;23:616-87.
18. Ellington MJ, Ganner M, Warner M, Cookson BD, Kearns AM. Polyclonal multiply antibiotic-resistant methicillin-resistant *Staphylococcus aureus* with Panton-Valentine leucocidin in England. *J Antimicrob Chemother* 2010;65:46-50.
19. Muder RR, Brennen C, Wagener MM, Vickers RM, Rihs JD, Hancock GA, et al. Methicillin-resistant staphylococcal colonization and infection in a long-term care facility. *Ann Intern Med* 1991;114:107-12.
20. Byrne FM, Wilcox MH. MRSA prevention strategies and current guidelines. *Injury* 2011;42 Suppl 5:S3-S6.
21. Gill SR, McIntyre LM, Nelson CL, Remortel B, Rude T, Reller LB, et al. Potential associations between severity of infection and the presence of virulence-associated genes in clinical strains of *Staphylococcus aureus*. *PLoS One* 2011;6:e18673.
22. Lyon BR, Skurray R. Antimicrobial resistance of *Staphylococcus aureus*: genetic basis. [Review] [559 refs]. *Microbiological Reviews* 1987;51:88-134.
23. Annear DI. The effect of temperature on resistance of *Staphylococcus aureus* to methicillin and some other antibiotics. *Med J Aust* 1968;1:444-6.
24. Hartman BJ, Tomasz A. Expression of methicillin resistance in heterogeneous strains of *Staphylococcus aureus*. *Antimicrob Agents Chemother* 1986;29:85-92.
25. Katayama Y, Ito T, Hiramatsu K. A new class of genetic element, staphylococcus cassette chromosome mec, encodes methicillin resistance in *Staphylococcus aureus*. *Antimicrob Agents Chemother* 2000;44:1549-55.
26. Ito T, Okuma K, Ma XX, Yuzawa H, Hiramatsu K. Insights on antibiotic resistance of *Staphylococcus aureus* from its whole genome: genomic island SCC. *Drug Resist Updat* 2003;6:41-52.
27. Ito T, Ma XX, Takeuchi F, Okuma K, Yuzawa H, Hiramatsu K. Novel type V staphylococcal cassette chromosome mec driven by a novel cassette chromosome recombinase, ccrC. *Antimicrob Agents Chemother* 2004;48:2637-51.
28. Garcia-Alvarez L, Holden MT, Lindsay H, Webb CR, Brown DF, Curran MD, et al. Methicillin-resistant *Staphylococcus aureus* with a novel mecA homologue in human and bovine populations in the UK and Denmark: a descriptive study. *Lancet Infect Dis* 2011;11:595-603.
29. Kluytmans-Vandenbergh MF, Kluytmans JA. Community-acquired methicillin-resistant *Staphylococcus aureus*: current perspectives. *Clin Microbiol Infect* 2006;12 Suppl 1:9-15.
30. Barg N, Chambers H, Kernodle D. Borderline susceptibility to antistaphylococcal penicillins is not conferred exclusively by the hyperproduction of beta-lactamase. *Antimicrob Agents Chemother* 1991;35:1975-9.

31. Thauvin-Eliopoulos C, Rice LB, Eliopoulos GM, Moellering RC, Jr. Efficacy of oxacillin and ampicillin-sulbactam combination in experimental endocarditis caused by beta-lactamase-hyperproducing *Staphylococcus aureus*. *Antimicrob Agents Chemother* 1990;34:728-32.
32. de Lencastre H, Tomasz A. Reassessment of the number of auxiliary genes essential for expression of high-level methicillin resistance in *Staphylococcus aureus*. *Antimicrob Agents Chemother* 1994;38:2590-8.
33. Appelbaum PC. The emergence of vancomycin-intermediate and vancomycin-resistant *Staphylococcus aureus*. *Clin Microbiol Infect* 2006;12 Suppl 1:16-23.
34. Gould IM, David MZ, Esposito S, Garau J, Lina G, Mazzei T, et al. New insights into methicillin-resistant *Staphylococcus aureus* (MRSA) pathogenesis, treatment and resistance. *Int J Antimicrob Agents* 2012;39:96-104.
35. Brown DF, Edwards Di, Hawkey PM, Morrison D, Ridgway GL, Towner KJ, et al. Guidelines for the laboratory diagnosis and susceptibility testing of methicillin-resistant *Staphylococcus aureus* (MRSA). *J Antimicrob Chemother* 2005;56:1000-18.
36. Wood W, Harvey G, Olson ES, Reid TM. Aztreonam selective agar for gram positive bacteria. *J Clin Pathol* 1993;46:769-71.
37. Morton CE, Holt HA. A problem encountered using staphylococcus/streptococcus supplement. *Med Lab Sci* 1989;46:72-3.
38. Jones EM, Bowker KE, Cooke R, Marshall RJ, Reeves DS, MacGowan AP. Salt tolerance of EMRSA-16 and its effect on the sensitivity of screening cultures. *J Hosp Infect* 1997;35:59-62.
39. Bruins MJ, Juffer P, Wolfhagen MJ, Ruijs GJ. Salt tolerance of methicillin-resistant and methicillin-susceptible *Staphylococcus aureus*. *J Clin Microbiol* 2007;45:682-3.
40. Morris K, Wilson C, Wilcox MH. Evaluation of chromogenic methicillin-resistant *Staphylococcus aureus* media: sensitivity versus turnaround time. *J Hosp Infect* 2012;81:20-4.
41. Yang HY, Suh JT, Lee HJ. Evaluation of commercial selective agars in screening for methicillin-resistant *Staphylococcus aureus*. *Ann Clin Lab Sci* 2010;40:252-6.
42. Danial J, Noel M, Templeton KE, Cameron F, Mathewson F, Smith M, et al. Real-time evaluation of an optimized real-time PCR assay versus Brilliance chromogenic MRSA agar for the detection of methicillin-resistant *Staphylococcus aureus* from clinical specimens. *J Med Microbiol* 2011;60:323-8.
43. Polisena J, Chen S, Cimon K, McGill S, Forward K, Gardam M. Clinical effectiveness of rapid tests for methicillin resistant *Staphylococcus aureus* (MRSA) in hospitalized patients: a systematic review. *BMC Infect Dis* 2011;11:336.
44. Roisin S, Laurent C, Nonhoff C, Deplano A, Hallin M, Byl B, et al. Positive predictive value of the Xpert MRSA assay diagnostic for universal patient screening at hospital admission: influence of the local ecology. *Eur J Clin Microbiol Infect Dis* 2012;31:873-80.
45. Akcam FZ, Tinaz GB, Kaya O, Tigli A, Ture E, Hosoglu S. Evaluation of methicillin resistance by cefoxitin disk diffusion and PBP2a latex agglutination test in mecA-positive *Staphylococcus aureus*, and comparison of mecA with femA, femB, femX positivities. *Microbiol Res* 2009;164:400-3.
46. Weist K, Cimbal AK, Lecke C, Kampf G, Ruden H, Vonberg RP. Evaluation of six agglutination tests for *Staphylococcus aureus* identification depending upon local prevalence of methicillin-resistant *S. aureus* (MRSA). *J Med Microbiol* 2006;55:283-90.

47. British Society for Antimicrobial Chemotherapy HISatICNA. Revised guidelines for the control of methicillin-resistant *Staphylococcus aureus* infection in hospitals. *J Hosp Infect* 1998;39:253-90.
48. Skov R, Smyth R, Clausen M, Larsen AR, Frimodt-Moller N, Olsson-Liljequist B, et al. Evaluation of a cefoxitin 30 microg disc on Iso-Sensitest agar for detection of methicillin-resistant *Staphylococcus aureus*. *J Antimicrob Chemother* 2003;52:204-7.
49. Cauwelier B, Gordts B, Descheemaecker P, Van Landuyt H. Evaluation of a disk diffusion method with cefoxitin (30 microg) for detection of methicillin-resistant *Staphylococcus aureus*. *Eur J Clin Microbiol Infect Dis* 2004;23:389-92.
50. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 9/99.
51. Department for transport. Transport of Infectious Substances, 2011 Revision 5. 2011.
52. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2013-2014. 2012.
53. Home Office. Anti-terrorism, Crime and Security Act. 2001 (as amended).
54. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive. 2013. p. 1-32
55. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office. 2003.
56. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive. 2005.
57. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive. 2008.
58. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. *MMWR Surveill Summ* 2012;61:1-102.
59. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002. 5th ed. HSE Books; 2002.
60. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books. 2002.
61. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books. 2002.
62. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books. 2003.
63. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets. 2000.
64. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 24-3-2005. p. 1-14

65. Baron EJ, Miller JM, Weinstein MP, Richter SS, Gilligan PH, Thomson RB, Jr., et al. A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2013 Recommendations by the Infectious Diseases Society of America (IDSA) and the American Society for Microbiology (ASM). *Clin Infect Dis* 2013;57:e22-e121.
66. Rishmawi N, Ghneim R, Kattan R, Ghneim R, Zoughbi M, Abu-Diab A, et al. Survival of fastidious and nonfastidious aerobic bacteria in three bacterial transport swab systems. *J Clin Microbiol* 2007;45:1278-83.
67. Barber S, Lawson PJ, Grove DI. Evaluation of bacteriological transport swabs. *Pathology* 1998;30:179-82.
68. Van Horn KG, Audette CD, Sebeck D, Tucker KA. Comparison of the Copan ESwab system with two Amies agar swab transport systems for maintenance of microorganism viability. *J Clin Microbiol* 2008;46:1655-8.
69. Nys S, Vijgen S, Magerman K, Cartuyvels R. Comparison of Copan eSwab with the Copan Venturi Transystem for the quantitative survival of *Escherichia coli*, *Streptococcus agalactiae* and *Candida albicans*. *Eur J Clin Microbiol Infect Dis* 2010;29:453-6.
70. Tano E, Melhus A. Evaluation of three swab transport systems for the maintenance of clinically important bacteria in simulated mono- and polymicrobial samples. *APMIS* 2011;119:198-203.
71. Smismans A, Verhaegen J, Schuermans A, Frans J. Evaluation of the Copan ESwab transport system for the detection of methicillin-resistant *Staphylococcus aureus*: a laboratory and clinical study. *Diagn Microbiol Infect Dis* 2009;65:108-11.
72. Human RP, Jones GA. Survival of bacteria in swab transport packs. *Med Lab Sci* 1986;43:14-8.
73. Widmer AF. Ceftobiprole: a new option for treatment of skin and soft-tissue infections due to methicillin-resistant *Staphylococcus aureus*. *Clin Infect Dis* 2008;46:656-8.
74. Public Health England. Laboratory Reporting to Public Health England: A Guide for Diagnostic Laboratories. 2013. p. 1-37.
75. Department of Health. Health Protection Legislation (England) Guidance. 2010. p. 1-112.
76. Scottish Government. Public Health (Scotland) Act. 2008 (as amended).
77. Scottish Government. Public Health etc. (Scotland) Act 2008. Implementation of Part 2: Notifiable Diseases, Organisms and Health Risk States. 2009.
78. The Welsh Assembly Government. Health Protection Legislation (Wales) Guidance. 2010.
79. Home Office. Public Health Act (Northern Ireland) 1967 Chapter 36. 1967 (as amended).