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
"NICE has renewed accreditation of the process used by **Public Health England (PHE)** to produce **UK Standards for Microbiology Investigations**. The renewed accreditation is valid until **30 June 2021** and applies to guidance produced using the processes described in **UK standards for microbiology investigations (UKSMIs) Development process, S9365', 2016**. The original accreditation term began in **July 2011**."

Emesso da Standards Unit, Microbiology Services, PHE

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Standards Unit
National Infection Service
Public Health England
61 Colindale Avenue
London NW9 5EQ

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UK Standards for Microbiology Investigations | Emesso da Standards Unit, Public Health England

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1. Cherry JD. Pertussis in young infants throughout the world. *Clinical Infectious Diseases* 2016;63:S119-S22. **A, III**

2. Zlomy M. Rediscovering Pertussis. *Front Pediatr* 2016;4:52. **A, III**

3. Cherry JD, Tan T, Wirsing Von Konig CH, Forsyth KD, Thisyakorn U, Greenberg D et al. Clinical definitions of pertussis: Summary of a global pertussis initiative roundtable meeting, february 2011. *Clinical Infectious Diseases* 2012;54:1756-64. **A, III**

4. Miyashita N, Akaike H, Teranishi H, Kawai Y, Ouchi K, Kato T et al. Diagnostic value of symptoms and laboratory data for pertussis in adolescent and adult patients. *BMC Infect Dis* 2013;13:129. **B, I**

5. Teepe J, Broekhuizen BD, Ieven M, Loens K, Huygen K, Kretzschmar M et al. Prevalence, diagnosis, and disease course of pertussis in adults with acute cough: a prospective, observational study in primary care.

British Journal of General Practice 2015;65:e662-7. **A, I**

6. Versteegh FGA. The re-emergence of pertussis. *Pediatric Pulmonology* 2014;49:S13-S4. **B, VI**

7. Public Health England. Pertussis: laboratory confirmed cases reported in England 2017. 2017.

8. Department of Health. Pregnant women to be offered whooping cough vaccination. 2012. **A, V**
9. Public Health England. Guidelines for the Public Health Management of Pertussis. 2016. **A, V**
10. Public Health England. Vaccination against pertussis (Whooping cough) for pregnant women. 2016. **A, V**
11. Nunes MC, Soofie N, Downs S, Tebeila N, Mudau A, De Gouveia L et al. Comparing the yield of nasopharyngeal swabs, nasal aspirates, and induced sputum for detection of bordetella pertussis in hospitalized infants. *Clinical Infectious Diseases* 2016;63:S181-S6. **A, II**
12. Curran T, Coyle PV. Understanding the true burden and infection dynamics of Bordetella pertussis using molecular diagnostics. *Journal of Infection* 2016;72:504-5. **B, VI**
13. Qin X, Zerr DM, Kronman MP, Adler AL, Berry JE, Rich S et al. Comparison of molecular detection methods for pertussis in children during a state-wide outbreak. *Annals of Clinical Microbiology and Antimicrobials* 2016;15 (1) (no pagination). **A, II**
14. Vaz-de-Lima LRA, Martin MD, Pawloski LC, Leite D, Rocha KCP, De Brito CA et al. Serodiagnosis as adjunct assay for pertussis infection in Sao Paulo, Brazil. *Clinical and Vaccine Immunology* 2014;21:636-40. **B, I**
15. Wirsing Von Konig CH. Pertussis diagnostics: Overview and impact of immunization. *Expert Review of Vaccines* 2014;13:1167-74. **B, III**
16. Bock JM, Burtis CC, Poetker DM, Blumin JH, Frank MO. Serum immunoglobulin G analysis to establish a delayed diagnosis of chronic cough due to Bordetella pertussis. *Otolaryngology - Head and Neck Surgery* 2012;146:63-7. **B, III**
17. DeVincenzo JP, Guyton C, Rea H, Elmore E, Patel S, Wynn L et al. Molecular detection and quantification of pertussis and correlation with clinical outcomes in children. *Diagnostic Microbiology and Infectious Disease* 2013;76:10-5. **B, II**
18. Qin X. Resurgence of Pertussis and Its Laboratory Diagnosis. *Clinical Microbiology Newsletter* 2015;37:69-76. **A, III**
19. Stone BL, Daly J, Srivastava R. Duration of Bordetella pertussis polymerase chain reaction positivity in confirmed pertussis illness. *Journal of the Pediatric Infectious Diseases Society* 2014;3:347-9. **A, I**
20. Zouari A, Smaoui H, Kechrid A. The diagnosis of pertussis: Which method to choose? *Critical Reviews in Microbiology* 2012;38:111-21. **A, III**
21. Bahtouee M, Moezzi SH, Nasirbagheban Z, Panahinia P. Evaluation of IgG serum level and PCR of respiratory secretion for Bordetella pertussis among adult suspected patients with whooping cough more than two weeks. *European Respiratory Journal Conference: European Respiratory Society Annual Congress* 2014;44. **B, II**
22. Baier M, Bohnert J, Hermann B, Pfister W. Comparison of three serological tests for pertussis antibody detection with respect to clinical data. *International Journal of Medical Microbiology* 2013;303:10. **B, VI**
23. Ruijs GJ, Groenendijk TW, Biever M. Shelf life of prepared Bordet-Gengou and Regan-Lowe agar plates for isolation of Bordetella pertussis. *EurJ Clin Microbiol Infect Dis* 1991;10:974-8. **B, II**

24. Hoppe JE, Schlagenhaut M. Comparison of three kinds of blood and two incubation atmospheres for cultivation of *Bordetella pertussis* on charcoal agar. *JClinMicrobiol* 1989;27:2115-7. **B, II**
25. Stauffer LR, Brown DR, Sandstrom RE. Cephalixin-supplemented Jones-Kendrick charcoal agar for selective isolation of *Bordetella pertussis*: comparison with previously described media. *JClinMicrobiol* 1983;17:60-2. **B, II**
26. Cloud JL, Hymas W, Carroll KC. Impact of nasopharyngeal swab types on detection of *Bordetella pertussis* by PCR and culture. *JClinMicrobiol* 2002;40:3838-40. **B, II**
27. European Parliament. UK Standards for Microbiology Investigations (UK SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes". 1998. **A, V**
28. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices 1998. 1-37. **A, V**
29. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 2009. **A, V**
30. Department for Transport. Transport of Infectious Substances, 2011 Revision 5. 2011. **A, V**
31. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2017-2018. 2017. **A, V**
32. Home Office. Anti-terrorism, Crime and Security Act. 2001. **A, V**
33. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive 2013. 1-35. **A, V**
34. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office 2003. **A, V**
35. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive 2005. **A, V**
36. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive 2008. **A, V**
37. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. *MMWR Surveill Summ* 2012;61:1-102. **B, IV**
38. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002 (as amended). HSE Books,. 2013. **A, V**
39. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books,. 2002. **A, V**

40. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books 2003. **A, V**
41. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets 2000. **A, V**
42. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 2005. 1-14. **A, V**
43. Department of Health. Transport of Infectious Substances. Best Practice Guidance for Microbiology Laboratories. Department of Health. 1-13. 2007. **A, V**
44. Baron EJ, Miller JM, Weinstein MP, Richter SS, Gilligan PH, Thomson RB, Jr. et al. A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2013 Recommendations by the Infectious Diseases Society of America (IDSA) and the American Society for Microbiology (ASM). *ClinInfectDis* 2013;57:e22-e121. **B, V**
45. Maiya S, Desai M, Baruah A, Weller P, Clarke JR, Gray J. Cough plate versus cough swab in patients with cystic fibrosis; a pilot study. *Arch Dis Child* 2004;89:577-9. **B, II**
46. Rishmawi N, Ghneim R, Kattan R, Ghneim R, Zoughbi M, Abu-Diab A et al. Survival of fastidious and nonfastidious aerobic bacteria in three bacterial transport swab systems. *JClinMicrobiol* 2007;45:1278-83. **B, II**
47. Barber S, Lawson PJ, Grove DI. Evaluation of bacteriological transport swabs. *Pathology* 1998;30:179-82. **C, II**
48. Van Horn KG, Audette CD, Sebeck D, Tucker KA. Comparison of the Copan ESwab system with two Amies agar swab transport systems for maintenance of microorganism viability. *JClinMicrobiol* 2008;46:1655-8. **B, II**
49. Nys S, Vijgen S, Magerman K, Cartuyvels R. Comparison of Copan eSwab with the Copan Venturi Transystem for the quantitative survival of *Escherichia coli*, *Streptococcus agalactiae* and *Candida albicans*. *EurJClinMicrobiolInfectDis* 2010;29:453-6. **B, II**
50. Tano E, Melhus A. Evaluation of three swab transport systems for the maintenance of clinically important bacteria in simulated mono- and polymicrobial samples. *APMIS* 2011;119:198-203. **B, II**
51. Crowcroft NS, Booy R, Harrison T, Spicer L, Britto J, Mok Q et al. Severe and unrecognised: pertussis in UK infants. *ArchDisChild* 2003;88:802-6. **B, II**
52. The Royal College of Pathologists. The retention and storage of pathological records and specimens (5th edition). 1-59. 2015. **A, V**
53. Bamberger ES, Srugo I. What is new in pertussis? *Eur J Pediatr* 2008;167:133-9. **B, II**
54. Cox HC, Jacob K, Whiley DM, Bletchly C, Nimmo GR, Nissen MD et al. Further evidence that the IS481 target is suitable for real-time PCR detection of *Bordetella pertussis*. *Pathology* 2013;45:202-3. **B, VI**
55. Gao F, Mahoney JC, Daly ER, Lamothe W, Tullo D, Bean C. Evaluation of a multitarget real-time PCR assay for detection of *Bordetella* species during a pertussis outbreak in New Hampshire in 2011. *Journal of clinical microbiology* 2014;52:302-6. **B, III**

56. Fry NK, Duncan J, Wagner K, Tzivra O, Doshi N, Litt DJ et al. Role of PCR in the diagnosis of pertussis infection in infants: 5 years' experience of provision of a same-day real-time PCR service in England and Wales from 2002 to 2007. *Journal of medical microbiology* 2009;58:1023-9. **B, I**
57. Loeffelholz M. Towards improved accuracy of *Bordetella pertussis* nucleic acid amplification tests. *Journal of clinical microbiology* 2012;50:2186-90. **A, II**
58. ECDC. Guidance and protocol for the use of real-time PCR in laboratory diagnosis of human infection with *Bordetella pertussis* or *Bordetella parapertussis*. 2012. **A, IV**
59. Tatti KM, Tondella ML. Utilization of multiple real-time PCR assays for the diagnosis of *Bordetella* spp. in clinical specimens. *Methods in Molecular Biology* 2013;943:135-47. **A, IV**
60. Fry NK, Litt DJ, Duncan J, Vaghji L, Warrener L, Samuel D et al. Modelling anti-pertussis toxin IgG antibody decay following primary and preschool vaccination with an acellular pertussis vaccine in UK subjects using a modified oral fluid assay. *Journal of medical microbiology* 2013;62:1281-9. **A, IV**
61. Tatti KM, Sparks KN, Boney KO, Tondella ML. Novel multitarget real-time PCR assay for rapid detection of *Bordetella* species in clinical specimens. *Journal of clinical microbiology* 2011;49:4059-66. **A, I**
62. Barbuddhe SB, Maier T, Schwarz G, Kostrzewa M, Hof H, Domann E et al. Rapid identification and typing of listeria species by matrix-assisted laser desorption ionization-time of flight mass spectrometry. *Appl Environ Microbiol* 2008;74:5402-7. **B, II**
63. ECDC. EQA scheme for *Bordetella* identification and *B.pertussis* typing. Stockholm: ECDC; 2014. **B, III**
64. Clark AE, Kaleta EJ, Arora A, Wolk DM. Matrix-assisted laser desorption ionization-time of flight mass spectrometry: a fundamental shift in the routine practice of clinical microbiology. *Clin Microbiol Rev* 2013;26:547-603. **B, II**
65. Xing D, Wirsing von Konig CH, Newland P, Riffelmann M, Meade BD, Corbel M et al. Characterization of reference materials for human antiserum to pertussis antigens by an international collaborative study. *Clin Vaccine Immunol* 2009;16:303-11. **A, I**
66. Lind-Brandberg L, Welinder-Olsson C, Lagergard T, Taranger J, Trollfors B, Zackrisson G. Evaluation of PCR for diagnosis of *Bordetella pertussis* and *Bordetella parapertussis* infections. *Journal of clinical microbiology* 1998;36:679-83. **A, I**
67. Public Health England. Laboratory Reporting to Public Health England: A Guide for Diagnostic Laboratories. Public Health England 2016. 1-29. **A, V**
68. Department of Health. Health Protection Legislation (England) Guidance. 1-112. 2010. **A, V**
69. Scottish Government. Public Health (Scotland) Act. 2008. **A, V**
70. Scottish Government. Public Health etc. (Scotland) Act 2008. Implementation of Part 2: Notifiable Diseases, Organisms and Health Risk States. 2009. **A, V**
71. The Welsh Assembly Government. Health Protection Legislation (Wales) Guidance. 2010. **A, V**
72. Home Office. Public Health Act (Northern Ireland) 1967 Chapter 36. 1967. **A, V**