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
Inviare a:

Staphylococcus Reference Laboratory Section
Laboratory of Healthcare-Associated Hospital Infection
Centre for Infections
Health Protection Agency
61 Colindale Avenue
London
NW9 5HT
<http://www.hpa.org.uk/cfi/lhcai/default.htm>

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IN REVISIONE

Bibliografia

1. Kloos W. Taxonomy and systematics of staphylococci indigenous to humans. In: Archer GL, Crossley KB, editors. *The Staphylococci in Human Disease*. New York: Churchill Livingstone; 1997. p. 113-7.
2. Kloos WE, Bannerman TL. Update on clinical significance of coagulase-negative staphylococci. *Clin Microbiol Rev* 1994;7:117-40.
3. Weinstein MP, Mirrett S, Van Pelt L, McKinnon M, Zimmer BL, Kloos W, et al. Clinical importance of identifying coagulase-negative staphylococci isolated from blood cultures: evaluation of MicroScan Rapid and Dried Overnight Gram-Positive panels versus a conventional reference method. *J Clin Microbiol* 1998;36:2089-92.
4. Holt JG. Gram-Positive cocci. In: Holt JG, Krieg NR, Sneath PHA, Staley JT, Williams ST, editors. *Bergey's Manual of Determinative Bacteriology*. 9th ed. Baltimore; 1994. p. 527-37.
5. Jarlov JO, Højbjerg T, Busch-Sørensen C, Scheibel J, Møller JK, Kolmos HJ, et al. Coagulase-negative Staphylococci in Danish blood cultures: species distribution and antibiotic susceptibility. *J Hosp Infect* 1996;32:217-27.
6. Vandenesch F, Lebeau C, Bes M, McDevitt D, Greenland T, Novick RP, et al. Coagulase deficiency in clinical isolates of *Staphylococcus aureus* involves both transcriptional and post-transcriptional defects. *J Med Microbiol* 1994;40:344-9.
7. Rupp ME, Archer GL. Coagulase-negative staphylococci: pathogens associated with medical progress. *Clin Infect Dis* 1994;19:231-43.
8. Christensen GD, Parisi JT, Bisno AL, Simpson WA, Beachey EH. Characterization of clinically significant strains of coagulase-negative staphylococci. *J Clin Microbiol* 1983;18:258-69.
9. Jansen B, Schumacher-Perdreau F, Peters G, Pulverer G. New aspects in the pathogenesis and prevention of polymer-associated foreign-body infections caused by coagulase-negative staphylococci. *J Invest Surg* 1989;2:361-80.
10. Bannerman TL, Wadiak DL, Kloos WE. Susceptibility of *Staphylococcus* species and subspecies to teicoplanin. *Antimicrob Agents Chemother* 1991;35:1919-22.
11. Vandenesch F, Perrier-Gros-Claude JD, Bes M, Fuhrmann C, Delorme V, Mouren C, et al. *Staphylococcus pasteuri*-specific oligonucleotide probes derived from a random amplified DNA fragment 259. *FEMS Microbiol Lett* 1995;132:147-52.
12. Baker JS. Comparison of various methods for differentiation of staphylococci and micrococci. *J Clin Microbiol* 1984;19:875-9.
13. Faller A, Schleifer KH. Modified oxidase and benzidine tests for separation of staphylococci from micrococci. *J Clin Microbiol* 1981;13:1031-5.
14. van Tiel FH, Slangen BF, Schouten HC, Jacobs JA. Study of *Stomatococcus mucilaginosus* isolated in a hospital ward using phenotypic characterization. *Eur J Clin Microbiol Infect Dis* 1995;14:193-8.
15. Personne P, Bes M, Lina G, Vandenesch F, Brun Y, Etienne J. Comparative performances of six agglutination kits assessed by using typical and atypical strains of *Staphylococcus aureus*. *J Clin Microbiol* 1997;35:1138-40.

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16. European Parliament. UK Standards for Microbiology Investigations (SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes".
17. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices. 7-12-1998, p. 1-37.
18. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 9/99.
19. Department for transport. Transport of Infectious Substances, 2011 Revision 5. 2011.
20. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2013-2014. 2012.
21. Home Office. Anti-terrorism, Crime and Security Act. 2001 (as amended).
22. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive. 2013. p. 1-32
23. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office. 2003.
24. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive. 2005.
25. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive. 2008.
26. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. *MMWR Surveill Summ* 2012;61:1-102.
27. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002. 5th ed. HSE Books; 2002.
28. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books. 2002.
29. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books. 2002.
30. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books. 2003.
31. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets. 2000.
32. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 24-3-2005. p. 1-14
33. Laughlin TJ, Armstrong DG, Caporusso J, Lavery LA. Soft tissue and bone infections from puncture wounds in children. *West J Med* 1997;166:126-8.

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34. Patel SR, Oleginski TP, Perruquet JL, Harrington TM. Pyomyositis: clinical features and predisposing conditions. *J Rheumatol* 1997;24:1734-8.
35. Santos KR, Fonseca LS, Bravo Neto GP, Gontijo Filho PP. Surgical site infection: rates, etiology and resistance patterns to antimicrobials among strains isolated at Rio de Janeiro University Hospital. *Infection* 1997;25:217-20.
36. Fayon MJ, Tucci M, Lacroix J, Farrell CA, Gauthier M, Lafleur L, et al. Nosocomial pneumonia and tracheitis in a pediatric intensive care unit: a prospective study. *Am J Respir Crit Care Med* 1997;155:162-9.
37. Abele-Horn M, Dauber A, Bauernfeind A, Russwurm W, Seyfarth-Metzger I, Gleich P, et al. Decrease in nosocomial pneumonia in ventilated patients by selective oropharyngeal decontamination (SOD). *Intensive Care Med* 1997;23:187-95.
38. Osterlund A, Nordlund E. Wound infection caused by *Staphylococcus hyicus* subspecies *hyicus* after a donkey bite. *Scand J Infect Dis* 1997;29:95.
39. Lee J. *Staphylococcus intermedius* isolated from dog-bite wounds. *J Infect* 1994;29:105.
40. Celard M, Vandenesch F, Darbas H, Grandjean J, Jean-Pierre H, Kirkorian G, et al. Pacemaker infection caused by *Staphylococcus schleiferi*, a member of the human preaxillary flora: four case reports. *Clin Infect Dis* 1997;24:1014-5.
41. Ozturkeri H, Kocabeyoglu O, Yergok YZ, Kosan E, Yenen OS, Keskin K. Distribution of coagulase-negative staphylococci, including the newly described species *Staphylococcus schleiferi*, in nosocomial and community acquired urinary tract infections. *Eur J Clin Microbiol Infect Dis* 1994;13:1076-9.
42. Latorre M, Rojo PM, Unzaga MJ, Cisterna R. *Staphylococcus schleiferi*: a new opportunistic pathogen. *Clin Infect Dis* 1993;16:589-90.
43. Mancao M, Miller C, Cochrane B, Hoff C, Sauter K, Weber E. Cerebrospinal fluid shunt infections in infants and children in Mobile, Alabama. *Acta Paediatr* 1998;87:667-70.
44. Ludlam H, Tremlett CH, Wilson AP. Preventing infection with *Staphylococcus aureus* in CAPD. *Perit Dial Int* 1997;17:405-6.
45. Henke PK, Bergamini TM, Garrison JR, Brittan KR, Peyton JC, Lam TM. *Staphylococcus epidermidis* graft infection is associated with locally suppressed major histocompatibility complex class II and elevated MAC-1 expression. *Arch Surg* 1997;132:894-902.
46. Karamanos NK, Syrokou A, Panagiotopoulou HS, Anastassiou ED, Dimitracopoulos G. The major 20-kDa polysaccharide of *Staphylococcus epidermidis* extracellular slime and its antibodies as powerful agents for detecting antibodies in blood serum and differentiating among slime-positive and -negative *S. epidermidis* and other staphylococci species. *Arch Biochem Biophys* 1997;342:389-95.
47. Schneider PF, Riley TV. *Staphylococcus saprophyticus* urinary tract infections: epidemiological data from Western Australia. *Eur J Epidemiol* 1996;12:51-4.
48. Mehta G, Kumari S. Multi-resistant *Staphylococcus haemolyticus* in a neonatal unit in New Delhi. *Ann Trop Paediatr* 1997;17:15-20.
49. Burnie JP, Naderi-Nasab M, Loudon KW, Matthews RC. An epidemiological study of blood culture isolates of coagulase-negative staphylococci demonstrating hospital-acquired infection. *J Clin Microbiol* 1997;35:1746-50.

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50. Akiyama H, Kanzaki H, Tada J, Arata J. Coagulase-negative staphylococci isolated from various skin lesions. *J Dermatol* 1998;25:563-8.
51. Kessler RB, Kimbrough RC, III, Jones SR. Infective endocarditis caused by *Staphylococcus hominis* after vasectomy. *Clin Infect Dis* 1998;27:216-7.
52. Barcs I, Herendi A, Lipcsey A, Bogнар C, Hashimoto H. Phage pattern and antibiotic resistance pattern of coagulase-negative staphylococci obtained from immunocompromised patients. *Microbiol Immunol* 1992;36:947-59.
53. Kloos WE, George CG, Olgiate JS, Van Pelt L, McKinnon ML, Zimmer BL, et al. *Staphylococcus hominis* subsp. *novobiosepticus* subsp. nov., a novel trehalose- and N-acetyl-D-glucosamine-negative, novobiocin- and multiple-antibiotic-resistant subspecies isolated from human blood cultures. *Int J Syst Bacteriol* 1998;48 Pt 3:799-812.
54. Jarlov JO, Prag J, Rosdahl VT, Espersen F. Evaluation of staphylococci isolated from a blood culture system (colorbact). *APMIS* 1995;103:383-7.
55. Crichton PB, Anderson LA, Phillips G, Davey PG, Rowley DI. Subspecies discrimination of staphylococci from revision arthroplasties by ribotyping. *J Hosp Infect* 1995;30:139-47.
56. al Rashdan A, Bashir R, Khan FA. *Staphylococcus capitis* causing aortic valve endocarditis. *J Heart Valve Dis* 1998;7:518-20.
57. Vandenesch F, Eykyn SJ, Bes M, Meugnier H, Fleurette J, Etienne J. Identification and ribotypes of *Staphylococcus caprae* isolates isolated as human pathogens and from goat milk. *J Clin Microbiol* 1995;33:888-92.
58. Shuttleworth R, Behme RJ, McNabb A, Colby WD. Human isolates of *Staphylococcus caprae*: association with bone and joint infections. *J Clin Microbiol* 1997;35:2537-41.
59. Schnitzler N, Meilicke R, Conrads G, Frank D, Haase G. *Staphylococcus lugdunensis*: report of a case of peritonitis and an easy-to-perform screening strategy. *J Clin Microbiol* 1998;36:812-3.
60. De Hondt G, Ieven M, Vanderersch C, Colaert J. Destructive endocarditis caused by *Staphylococcus lugdunensis*. Case report and review of the literature. *Acta Clin Belg* 1997;52:27-30.
61. Koh TW, Brecker SJ, Layton CA. Successful treatment of *Staphylococcus lugdunensis* endocarditis complicated by multiple emboli: a case report and review of the literature. *Int J Cardiol* 1996;55:193-7.
62. Gidley PW, Ghorayeb BY, Stiernberg CM. Contemporary management of deep neck space infections. *Otolaryngol Head Neck Surg* 1997;116:16-22.
63. Orrett FA, Shurland SM. Significance of coagulase-negative staphylococci in urinary tract infections in a developing country. *Conn Med* 1998;62:199-203.
64. Kolawole DO, Shittu AO. Unusual recovery of animal staphylococci from septic wounds of hospital patients in Ile-Ife, Nigeria. *Lett Appl Microbiol* 1997;24:87-90.
65. Buttery JP, Easton M, Pearson SR, Hogg GG. Pediatric bacteremia due to *Staphylococcus warneri*: microbiological, epidemiological, and clinical features. *J Clin Microbiol* 1997;35:2174-7.
66. Peces R, Gago E, Tejada F, Lares AS, Alvarez-Grande J. Relapsing bacteraemia due to *Micrococcus luteus* in a haemodialysis patient with a Perm-Cath catheter. *Nephrol Dial Transplant* 1997;12:2428-9.

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67. Rabitsch W, Brugger SA, Pirker W, Baumgartner C, Reiter E, Keil F, et al. Symmetrical necrosis of globus pallidus with severe gait disturbance in a patient with myelodysplastic syndrome given allogeneic marrow transplantation. *Ann Hematol* 1997;75:235-7.
68. Kern W, Kurrle E, Vanek E. Ofloxacin for prevention of bacterial infections in granulocytopenic patients. *Infection* 1987;15:427-33.
69. Kiehn TE, Armstrong D. Changes in the spectrum of organisms causing bacteremia and fungemia in immunocompromised patients due to venous access devices. *Eur J Clin Microbiol Infect Dis* 1990;9:869-72.
70. Abraham J, Bilgrami S, Dorsky D, Edwards RL, Feingold J, Hill DR, et al. *Stomatococcus mucilaginosus* meningitis in a patient with multiple myeloma following autologous stem cell transplantation. *Bone Marrow Transplant* 1997;19:639-41.
71. Park MK, Khan J, Stock F, Lucey DR. Successful treatment of *Stomatococcus mucilaginosus* meningitis with intravenous vancomycin and intravenous ceftriaxone. *Clin Infect Dis* 1997;24:278.
72. McWhinney PH, Kibbler CC, Gillespie SH, Patel S, Morrison D, Hoffbrand AV, et al. *Stomatococcus mucilaginosus*: an emerging pathogen in neutropenic patients. *Clin Infect Dis* 1992;14:641-6.
73. Vasishtha S, Isenberg HD, Sood SK. *Gemella morbillorum* as a cause of septic shock. *Clin Infect Dis* 1996;22:1084-6.
74. Kloos W, Bannerman TL. *Staphylococcus* and *Micrococcus*. In: Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover RH, editors. *Manual of Clinical Microbiology*. 7th ed. Washington DC: American Society for Microbiology; 1999. p. 264-82.
75. Public Health England. *Laboratory Reporting to Public Health England: A Guide for Diagnostic Laboratories*. 2013. p. 1-37.
76. Department of Health. *Health Protection Legislation (England) Guidance*. 2010. p. 1-112.
77. Scottish Government. *Public Health (Scotland) Act*. 2008 (as amended).
78. Scottish Government. *Public Health etc. (Scotland) Act 2008. Implementation of Part 2: Notifiable Diseases, Organisms and Health Risk States*. 2009.
79. The Welsh Assembly Government. *Health Protection Legislation (Wales) Guidance*. 2010.
80. Home Office. *Public Health Act (Northern Ireland) 1967 Chapter 36*. 1967 (as amended).

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16. European Parliament. UK Standards for Microbiology Investigations (SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes".
17. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices. 7-12-1998, p. 1-37.
18. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 9/99.
19. Department for transport. Transport of Infectious Substances, 2011 Revision 5. 2011.
20. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2013-2014. 2012.
21. Home Office. Anti-terrorism, Crime and Security Act. 2001 (as amended).
22. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive. 2013. p. 1-32
23. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office. 2003.
24. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive. 2005.
25. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive. 2008.
26. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. MMWR Surveill Summ 2012;61:1-102.
27. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002. 5th ed. HSE Books; 2002.
28. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books. 2002.
29. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books. 2002.
30. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books. 2003.
31. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets. 2000.
32. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 24-3-2005. p. 1-14
33. Laughlin TJ, Armstrong DG, Caporusso J, Lavery LA. Soft tissue and bone infections from puncture wounds in children. West J Med 1997;166:126-8.

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34. Patel SR, Oleginski TP, Perruquet JL, Harrington TM. Pyomyositis: clinical features and predisposing conditions. *J Rheumatol* 1997;24:1734-8.
35. Santos KR, Fonseca LS, Bravo Neto GP, Gontijo Filho PP. Surgical site infection: rates, etiology and resistance patterns to antimicrobials among strains isolated at Rio de Janeiro University Hospital. *Infection* 1997;25:217-20.
36. Fayon MJ, Tucci M, Lacroix J, Farrell CA, Gauthier M, Lafleur L, et al. Nosocomial pneumonia and tracheitis in a pediatric intensive care unit: a prospective study. *Am J Respir Crit Care Med* 1997;155:162-9.
37. Abele-Horn M, Dauber A, Bauernfeind A, Russwurm W, Seyfarth-Metzger I, Gleich P, et al. Decrease in nosocomial pneumonia in ventilated patients by selective oropharyngeal decontamination (SOD). *Intensive Care Med* 1997;23:187-95.
38. Osterlund A, Nordlund E. Wound infection caused by *Staphylococcus hyicus* subspecies *hyicus* after a donkey bite. *Scand J Infect Dis* 1997;29:95.
39. Lee J. *Staphylococcus intermedius* isolated from dog-bite wounds. *J Infect* 1994;29:105.
40. Celard M, Vandenesch F, Darbas H, Grando J, Jean-Pierre H, Kirkorian G, et al. Pacemaker infection caused by *Staphylococcus schleiferi*, a member of the human preaxillary flora: four case reports. *Clin Infect Dis* 1997;24:1014-5.
41. Ozturkeri H, Kocabeyoglu O, Yergok YZ, Kosan E, Yenen OS, Keskin K. Distribution of coagulase-negative staphylococci, including the newly described species *Staphylococcus schleiferi*, in nosocomial and community acquired urinary tract infections. *Eur J Clin Microbiol Infect Dis* 1994;13:1076-9.
42. Latorre M, Rojo PM, Unzaga MJ, Cisterna R. *Staphylococcus schleiferi*: a new opportunistic pathogen. *Clin Infect Dis* 1993;16:589-90.
43. Mancao M, Miller C, Cochrane B, Hoff C, Sauter K, Weber E. Cerebrospinal fluid shunt infections in infants and children in Mobile, Alabama. *Acta Paediatr* 1998;87:667-70.
44. Ludlam H, Tremlett CH, Wilson AP. Preventing infection with *Staphylococcus aureus* in CAPD. *Perit Dial Int* 1997;17:405-6.
45. Henke PK, Bergamini TM, Garrison JR, Brittan KR, Peyton JC, Lam TM. *Staphylococcus epidermidis* graft infection is associated with locally suppressed major histocompatibility complex class II and elevated MAC-1 expression. *Arch Surg* 1997;132:894-902.
46. Karamanos NK, Syrokou A, Panagiotopoulou HS, Anastassiou ED, Dimitracopoulos G. The major 20-kDa polysaccharide of *Staphylococcus epidermidis* extracellular slime and its antibodies as powerful agents for detecting antibodies in blood serum and differentiating among slime-positive and -negative *S. epidermidis* and other staphylococci species. *Arch Biochem Biophys* 1997;342:389-95.
47. Schneider PF, Riley TV. *Staphylococcus saprophyticus* urinary tract infections: epidemiological data from Western Australia. *Eur J Epidemiol* 1996;12:51-4.
48. Mehta G, Kumari S. Multi-resistant *Staphylococcus haemolyticus* in a neonatal unit in New Delhi. *Ann Trop Paediatr* 1997;17:15-20.
49. Burnie JP, Naderi-Nasab M, Loudon KW, Matthews RC. An epidemiological study of blood culture isolates of coagulase-negative staphylococci demonstrating hospital-acquired infection. *J Clin Microbiol* 1997;35:1746-50.

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